

## Murine cardiac troponin I peptide TFA

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Murine cardiac troponin I peptide (Murine cardiac TnI peptide) TFA is a peptide derived from murine cardiac troponin I (cTnI). It can trigger severe myocardial inflammation in mouse models, leading to fibrosis, heart failure, and increased mortality. This compound is utilized in heart failure research.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481