

## trans-2-Nonenoyl-CoA

## Chemical Properties

CAS No. :	1527491-63-9
Formula:	C30H50N7O17P3S
Molecular Weight:	905.74
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	trans-2-Nonenoyl-CoA (trans-2-Nonenoyl-coenzyme A) is a derivative of coenzyme A.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1041 mL	5.5203 mL	11.0407 mL
5 mM	0.2208 mL	1.1041 mL	2.2081 mL
10 mM	0.1104 mL	0.552 mL	1.1041 mL
50 mM	0.0221 mL	0.1104 mL	0.2208 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481