

dsP53-285 saRNA

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	dsP53-285 saRNA is a small activating RNA (saRNA) that can effectively stimulate the expression of wild-type p53 by targeting its promoter. This activation leads to the inhibition of growth and metastasis in bladder cancer cells.
Targets(IC50)	MDM-2/p53

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481