

Gly-Gly-Gly-PEG3-TCO TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Gly-Gly-Gly-PEG3-TCO TFA is a cleavable ADC linker incorporating a three-unit PEG structure, used in the synthesis of antibody-drug conjugates (ADC).
Targets(IC50)	ADC Linker
In vitro	An ADC (antibody-drug conjugate) is formed by connecting an antibody with a cytotoxic ADC drug through an ADC linker.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481