

[S5K, F6Y, L9mL, M10Mox] neurokinin A (4-10)

Chemical Properties

CAS No. : 2770688-17-8
Formula: C38H63N9O11
Molecular Weight: 821.97
Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	[S5K, F6Y, L9mL, M10Mox] neurokinin A (4-10) is a neuropeptide A analogue and peptide fragment. It exhibits high selectivity for the NK1R on both mouse and human tachykinin receptors, with a half-life of 10.3 hours in mice.
Targets(IC50)	Neurokinin receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2166 mL	6.0829 mL	12.1659 mL
5 mM	0.2433 mL	1.2166 mL	2.4332 mL
10 mM	0.1217 mL	0.6083 mL	1.2166 mL
50 mM	0.0243 mL	0.1217 mL	0.2433 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481