

HWGMWSY

Chemical Properties

CAS No. :	348082-59-7
Formula:	C47H55N11O10S
Molecular Weight:	966.08
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	HWGMWSY is a peptide associated with polystyrene. Although it is not targeted to any specific molecule, HWGMWSY interacts with the surface of polystyrene plastic.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.0351 mL	5.1756 mL	10.3511 mL
5 mM	0.207 mL	1.0351 mL	2.0702 mL
10 mM	0.1035 mL	0.5176 mL	1.0351 mL
50 mM	0.0207 mL	0.1035 mL	0.207 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481