

## (-)-Mitorubrinol

## Chemical Properties

CAS No. :	3215-47-2
Formula:	C <sub>21</sub> H <sub>18</sub> O <sub>8</sub>
Molecular Weight:	398.37
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	(-)-Mitorubrinol is found in the marine-derived strain <i>Penicillium purpurogenum</i> G59, which was mutated using diethyl sulfate (DES). It exhibits inhibitory effects on cancer cell lines K562, HL-60, HeLa, and BGC-823.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5102 mL	12.5511 mL	25.1023 mL
5 mM	0.502 mL	2.5102 mL	5.0205 mL
10 mM	0.251 mL	1.2551 mL	2.5102 mL
50 mM	0.0502 mL	0.251 mL	0.502 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481