

Apixaban impurity 18

Chemical Properties

CAS No. :	1643330-62-4
Formula:	C ₂₀ H ₂₆ ClN ₃ O ₃
Molecular Weight:	391.9
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Apixaban impurity 18 is an impurity of Apixaban.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5517 mL	12.7584 mL	25.5167 mL
5 mM	0.5103 mL	2.5517 mL	5.1033 mL
10 mM	0.2552 mL	1.2758 mL	2.5517 mL
50 mM	0.051 mL	0.2552 mL	0.5103 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481