

Ticagrelor impurity 5

Chemical Properties

CAS No. :	1863036-49-0
Formula:	C ₂₆ H ₃₂ F ₂ N ₆ O ₅ S
Molecular Weight:	578.64
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Ticagrelor impurity 5 is an impurity of Ticagrelor.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7282 mL	8.641 mL	17.2819 mL
5 mM	0.3456 mL	1.7282 mL	3.4564 mL
10 mM	0.1728 mL	0.8641 mL	1.7282 mL
50 mM	0.0346 mL	0.1728 mL	0.3456 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481