


Lornoxicam impurity 16

Chemical Properties

CAS No. :	70374-37-7
Formula:	C7H8ClNO4S2
Molecular Weight:	269.71
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lornoxicam Impurity 16 is an impurity of Lornoxicam.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.7077 mL	18.5384 mL	37.0769 mL
5 mM	0.7415 mL	3.7077 mL	7.4154 mL
10 mM	0.3708 mL	1.8538 mL	3.7077 mL
50 mM	0.0742 mL	0.3708 mL	0.7415 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481