

Neratinib impurity13

Chemical Properties

CAS No. :	1144516-13-1
Formula:	C ₂₆ H ₂₀ ClN ₅ O ₅
Molecular Weight:	517.93
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Neratinib impurity 13 is an impurity of Neratinib.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9308 mL	9.6538 mL	19.3076 mL
5 mM	0.3862 mL	1.9308 mL	3.8615 mL
10 mM	0.1931 mL	0.9654 mL	1.9308 mL
50 mM	0.0386 mL	0.1931 mL	0.3862 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481