

Valenclin impurity 1

Chemical Properties

CAS No. :	230615-53-9
Formula:	C ₁₃ H ₁₁ F ₃ N ₂ O ₃
Molecular Weight:	300.24
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Valenclin impurity 1 is a type of Valenclin impurity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3307 mL	16.6533 mL	33.3067 mL
5 mM	0.6661 mL	3.3307 mL	6.6613 mL
10 mM	0.3331 mL	1.6653 mL	3.3307 mL
50 mM	0.0666 mL	0.3331 mL	0.6661 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481