

## Abiraterone impurity 1

## Chemical Properties

CAS No. :	165334-72-5
Formula:	C <sub>24</sub> H <sub>31</sub> NO
Molecular Weight:	349.52
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Abiraterone impurity 1 is an impurity of Abiraterone.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8611 mL	14.3053 mL	28.6107 mL
5 mM	0.5722 mL	2.8611 mL	5.7221 mL
10 mM	0.2861 mL	1.4305 mL	2.8611 mL
50 mM	0.0572 mL	0.2861 mL	0.5722 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481