

ALN-3133 sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ALN-3133 sodium is a siRNA targeting VEGF and is one of the siRNAs in ALN-VSP. ALN-VSP is a lipid nanoparticle formulation containing two siRNAs against kinesin spindle protein (KSP) and VEGF, offering potential antitumor activity.
Targets(IC50)	VEGFR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481