

Apixaban impurity 2

Chemical Properties

CAS No. :	2187409-01-2
Formula:	C ₂₅ H ₂₈ N ₆ O ₄
Molecular Weight:	476.54
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Apixaban impurity 2 is an impurity of Apixaban.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0985 mL	10.4923 mL	20.9846 mL
5 mM	0.4197 mL	2.0985 mL	4.1969 mL
10 mM	0.2098 mL	1.0492 mL	2.0985 mL
50 mM	0.042 mL	0.2098 mL	0.4197 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481