

## Dabigatran impurity 1

## Chemical Properties

CAS No. :	2417628-79-4
Formula:	C <sub>25</sub> H <sub>24</sub> N <sub>6</sub> O <sub>4</sub>
Molecular Weight:	472.51
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Dabigatran impurity 1 (Dabigatran Carboxamide) is an impurity of Dabigatran.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1164 mL	10.5818 mL	21.1636 mL
5 mM	0.4233 mL	2.1164 mL	4.2327 mL
10 mM	0.2116 mL	1.0582 mL	2.1164 mL
50 mM	0.0423 mL	0.2116 mL	0.4233 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481