

Adrenaline impurity 3

Chemical Properties

CAS No. :	317351-39-6
Formula:	C16H18ClNO3
Molecular Weight:	307.77
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Adrenaline impurity 3 is a known impurity of Adrenaline.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2492 mL	16.2459 mL	32.4918 mL
5 mM	0.6498 mL	3.2492 mL	6.4984 mL
10 mM	0.3249 mL	1.6246 mL	3.2492 mL
50 mM	0.065 mL	0.3249 mL	0.6498 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481