

Palbociclib impurity 20

Chemical Properties

CAS No. :	1774898-55-3
Formula:	C ₁₄ H ₂₀ N ₄ O ₄
Molecular Weight:	308.34
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Palbociclib impurity 20 is an impurity of Palbociclib.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2432 mL	16.2159 mL	32.4317 mL
5 mM	0.6486 mL	3.2432 mL	6.4863 mL
10 mM	0.3243 mL	1.6216 mL	3.2432 mL
50 mM	0.0649 mL	0.3243 mL	0.6486 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481