

Apixaban impurity 1

Chemical Properties

CAS No. :	27143-07-3
Formula:	C ₁₁ H ₁₃ ClN ₂ O ₃
Molecular Weight:	256.69
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Apixaban impurity 1 is a byproduct found in Apixaban.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8957 mL	19.4787 mL	38.9575 mL
5 mM	0.7791 mL	3.8957 mL	7.7915 mL
10 mM	0.3896 mL	1.9479 mL	3.8957 mL
50 mM	0.0779 mL	0.3896 mL	0.7791 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481