


Empagliflozin impurity 83

Chemical Properties

CAS No. :	915095-84-0
Formula:	C ₁₇ H ₁₄ BrClO ₃
Molecular Weight:	381.65
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Empagliflozin impurity 83 is a variant of the compound Empagliflozin.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6202 mL	13.101 mL	26.202 mL
5 mM	0.524 mL	2.6202 mL	5.2404 mL
10 mM	0.262 mL	1.3101 mL	2.6202 mL
50 mM	0.0524 mL	0.262 mL	0.524 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481