

N-Acetyl-D-allo-isoleucine

Chemical Properties

CAS No. :	54831-20-8
Formula:	C ₈ H ₁₅ NO ₃
Molecular Weight:	173.21
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	N-Acetyl-D-allo-isoleucine is a derivative of amino acids and holds potential for research in amino acid metabolism.
Targets(IC50)	Amino Acids and Derivatives

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.7733 mL	28.8667 mL	57.7334 mL
5 mM	1.1547 mL	5.7733 mL	11.5467 mL
10 mM	0.5773 mL	2.8867 mL	5.7733 mL
50 mM	0.1155 mL	0.5773 mL	1.1547 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481