

AGN-745

## Chemical Properties

CAS No. : 1255908-66-7

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	AGN-745 is a chemically modified siRNA designed to target vascular endothelial growth factor receptor-1 (VEGFR-1), and it is utilized in the research of macular degeneration.
Targets(IC50)	VEGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481