

## 3-Chlorophenylacetylhydrazide

## Chemical Properties

CAS No. :	66464-86-6
Formula:	C <sub>8</sub> H <sub>9</sub> ClN <sub>2</sub> O
Molecular Weight:	184.623
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	3-Chlorophenylacetylhydrazide is a biochemical reagent that serves as a biomaterial or organic compound useful in life science research.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.4165 mL	27.0827 mL	54.1653 mL
5 mM	1.0833 mL	5.4165 mL	10.8331 mL
10 mM	0.5417 mL	2.7083 mL	5.4165 mL
50 mM	0.1083 mL	0.5417 mL	1.0833 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481