

Abiraterone impurity 15

Chemical Properties

CAS No. :	2141996-19-0
Formula:	C ₂₆ H ₃₅ NO ₂
Molecular Weight:	393.57
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Abiraterone Impurity 15 is a byproduct associated with Abiraterone.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5408 mL	12.7042 mL	25.4084 mL
5 mM	0.5082 mL	2.5408 mL	5.0817 mL
10 mM	0.2541 mL	1.2704 mL	2.5408 mL
50 mM	0.0508 mL	0.2541 mL	0.5082 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481