

Afovirsen

Chemical Properties

CAS No. : 151356-08-0

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Afovirsen is a phosphorothioate oligonucleotide. It is complementary to the mRNA sequence of the initiation codon for the E2 protein, which is essential for the replication of HPV types 6 and 11.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481