


GABAA receptor modulator-12

Chemical Properties

CAS No. :	6252-98-8
Formula:	C ₁₁ H ₁₀ N ₂ O ₂
Molecular Weight:	202.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	GABAA receptor modulator-12 (Compound 5) is a modulator of the GABAA receptor. It inhibits anesthetic-induced desensitization of the GABAA α 3 β 3 γ 2 receptor, with an IC ₅₀ of 0.2 μ M. This compound is useful for research in central nervous system diseases.
Targets(IC ₅₀)	GABA Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.9454 mL	24.7268 mL	49.4535 mL
5 mM	0.9891 mL	4.9454 mL	9.8907 mL
10 mM	0.4945 mL	2.4727 mL	4.9454 mL
50 mM	0.0989 mL	0.4945 mL	0.9891 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481