

Carbomer 676

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Carbomer 676 is an acrylic polymer that acts as a thickening agent. It forms high-viscosity gels, offering a stable gel matrix and demonstrating excellent resistance to heat, light, and microbial contamination. Carbomer 676 facilitates drug release and topical application, and is utilized in the pharmaceutical industry. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481