

Dabigatran impurity 80

Chemical Properties

CAS No. :	1416446-40-6
Formula:	C30H33N7O5
Molecular Weight:	571.64
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Dabigatran impurity 80 is an impurity of Dabigatran.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7494 mL	8.7468 mL	17.4935 mL
5 mM	0.3499 mL	1.7494 mL	3.4987 mL
10 mM	0.1749 mL	0.8747 mL	1.7494 mL
50 mM	0.035 mL	0.1749 mL	0.3499 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481