

Obudanersen sodium

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Obudanersen sodium is an antisense oligonucleotide designed to target and inhibit the expression of UBE3A antisense transcript (UBE3A-AS), and it is used in research on Angelman syndrome.
Targets(IC50)	E1/E2/E3 Enzyme

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481