

## Afatinib impurity 92

## Chemical Properties

CAS No. :	2413212-07-2
Formula:	C18H17FN4O2
Molecular Weight:	340.36
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Afatinibimpurity 92 is a type of Afatinib impurity.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9381 mL	14.6903 mL	29.3807 mL
5 mM	0.5876 mL	2.9381 mL	5.8761 mL
10 mM	0.2938 mL	1.469 mL	2.9381 mL
50 mM	0.0588 mL	0.2938 mL	0.5876 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481