

3-Chloropropanoic acid 6-(4-hydroxyphenoxy)hexyl ester

Chemical Properties

CAS No. :	2087904-32-1
Formula:	C ₁₅ H ₂₁ ClO ₄
Molecular Weight:	300.78
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	3-Chloropropanoic acid 6-(4-hydroxyphenoxy)hexyl ester is a pharmaceutical impurity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3247 mL	16.6234 mL	33.2469 mL
5 mM	0.6649 mL	3.3247 mL	6.6494 mL
10 mM	0.3325 mL	1.6623 mL	3.3247 mL
50 mM	0.0665 mL	0.3325 mL	0.6649 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481