

Arolol impurity 4

Chemical Properties

CAS No. :	68257-89-6
Formula:	C7H7NO2S
Molecular Weight:	169.2
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Arolol impurity 4 is a pharmaceutical impurity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.9102 mL	29.5508 mL	59.1017 mL
5 mM	1.182 mL	5.9102 mL	11.8203 mL
10 mM	0.591 mL	2.9551 mL	5.9102 mL
50 mM	0.1182 mL	0.591 mL	1.182 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481