

Boc-PEG6-Val-Cit-PAB-OH

Chemical Properties

CAS No. :	2055024-53-6
Formula:	C ₃₈ H ₆₆ N ₆ O ₁₃
Molecular Weight:	814.98
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Boc-PEG6-Val-Cit-PAB-OH is an ADC linker utilized in the synthesis of antibody-drug conjugates (ADC). It is applicable in cancer research.
Targets(IC50)	ADC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.227 mL	6.1351 mL	12.2702 mL
5 mM	0.2454 mL	1.227 mL	2.454 mL
10 mM	0.1227 mL	0.6135 mL	1.227 mL
50 mM	0.0245 mL	0.1227 mL	0.2454 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481