

## Tetrabenazine impurity 1

## Chemical Properties

CAS No. :	99672-64-7
Formula:	C19H27NO3
Molecular Weight:	317.43
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Tetrabenazine impurity 1 is an impurity of Tetrabenazine.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1503 mL	15.7515 mL	31.503 mL
5 mM	0.6301 mL	3.1503 mL	6.3006 mL
10 mM	0.315 mL	1.5752 mL	3.1503 mL
50 mM	0.063 mL	0.315 mL	0.6301 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481