

Apazunersen

Chemical Properties

CAS No. : 2919324-01-7

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Apazunersen is an antisense oligonucleotide designed to target and inhibit the expression of UBE3A antisense transcript (UBE3A-AS). This inhibition prevents the silencing of paternal alleles of the UBE3A gene, thereby reactivating the expression of the defective protein. It is utilized in the research of Angelman syndrome.
Targets(IC50)	E1/E2/E3 Enzyme

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481