

Sapablursen sodium

Chemical Properties

CAS No. : 2273008-00-5

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sapablursen sodium is an antisense oligonucleotide that reduces the production of TMPRSS6, leading to increased expression of serum hepcidin, which is a key regulator of iron homeostasis. Sapablursen may be used in research on blood disorders such as polycythemia vera.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481