

## Bovine CB12 of type II collagen

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Bovine CB12 of type II collagen (Bovine cyanogen bromide digested peptide 12) is a fragment of cyanogen bromide-digested type II collagen peptide (CB-peptide) useful for investigating B cell and T cell epitopes of type II collagen. This product is free from pepsin, making it ideal for studying T cell epitope specificity through proliferation assays. Additionally, Bovine CB12 of type II collagen serves as an antigen for evaluating antibody epitope specificity via ELISA analysis.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481