

Echiumiline

Chemical Properties

CAS No. :	174285-73-5
Formula:	C ₂₀ H ₃₁ N ₀₇
Molecular Weight:	397.47
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Echiumiline is an alkaloid that targets liver cytochrome P450 enzymes (CYP450). It induces DNA cross-linking and oxidative stress in liver cells, which leads to liver necrosis and fibrosis. Echiumiline holds potential for research in liver diseases.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5159 mL	12.5796 mL	25.1591 mL
5 mM	0.5032 mL	2.5159 mL	5.0318 mL
10 mM	0.2516 mL	1.258 mL	2.5159 mL
50 mM	0.0503 mL	0.2516 mL	0.5032 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481