

Clopidogrel impurity 3

Chemical Properties

CAS No. :	1015247-87-6
Formula:	C ₁₄ H ₁₄ ClNO ₂ S
Molecular Weight:	295.78
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Clopidogrelimpurity 3 is an impurity of Clopidogrel.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3809 mL	16.9045 mL	33.8089 mL
5 mM	0.6762 mL	3.3809 mL	6.7618 mL
10 mM	0.3381 mL	1.6904 mL	3.3809 mL
50 mM	0.0676 mL	0.3381 mL	0.6762 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481