

## Umedaptanib pegol

### Chemical Properties

CAS No. : 2734771-59-4

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Umedaptanib pegol is a nucleic acid aptamer that targets and inhibits fibroblast growth factor-2 (FGF2). It is currently under investigation for its potential use in the treatment of neovascular age-related macular degeneration (nAMD).
Targets(IC50)	FGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481