

## Chick CB12 of type II collagen

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Chick CB12 of type II collagen (Chick cyanogen bromide digested peptide 12) is a type II collagen cyanogen bromide digested peptide fragment (CB-peptide) utilized for studying B cell and T cell epitopes of type II collagen. This product is pepsin-free, making it ideal for investigating T cell epitope specificity through proliferation assays. Additionally, Chick CB12 of type II collagen serves as an antigen for analyzing antibody epitope specificity using ELISA.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481