

Obudanersen

Chemical Properties

CAS No. : 2834724-31-9

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Obudanersen is an antisense oligonucleotide designed to target and inhibit the expression of the UBE3A antisense transcript (UBE3A-AS) and is used in the study of Angelman syndrome.
Targets(IC50)	E1/E2/E3 Enzyme

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481