

## Guinea Pig serum albumin

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Guinea Pig Serum Albumin is a plasma protein derived from guinea pigs. This multifunctional protein exhibits an exceptional capacity for binding ligands, enabling it to serve as a transport molecule for a variety of metabolites, drugs, nutrients, metals, and other molecules.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481