

Obeversen sodium

Chemical Properties

CAS No. : 2304711-35-9

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Obeversen sodium is an antisense oligonucleotide that inhibits the synthesis of diacylglycerol acyltransferase 2 (DGAT-2). It is used in research related to non-alcoholic fatty liver disease.
Targets(IC50)	Acyltransferase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481