

## Pegpleranib

### Chemical Properties

CAS No. : 1618657-13-8

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Pegpleranib is a polyethylene glycol-modified nucleic acid aptamer that specifically binds to both homodimers and heterodimers of PDGF-BB and PDGF-AB, thereby inhibiting their interaction with associated tyrosine kinase receptors.
Targets(IC50)	PDGFR

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481