

Edifoligide

Chemical Properties

CAS No. : 328538-04-1

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Edifoligide is a decoy oligonucleotide that binds to the E2F transcription factor, inhibiting its activity, which may help prevent neointimal hyperplasia and vein graft failure.
Targets(IC50)	Early 2 Factor (E2F)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481