

## Dacomitinib impurity 12

## Chemical Properties

CAS No. :	2803451-99-0
Formula:	C <sub>15</sub> H <sub>12</sub> Cl <sub>2</sub> N <sub>4</sub> O
Molecular Weight:	335.19
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Dacomitinib impurity 12 is an impurity of Dacomitinib.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9834 mL	14.9169 mL	29.8338 mL
5 mM	0.5967 mL	2.9834 mL	5.9668 mL
10 mM	0.2983 mL	1.4917 mL	2.9834 mL
50 mM	0.0597 mL	0.2983 mL	0.5967 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481