

Rebamipide methyl ester

Chemical Properties

CAS No. :	90098-38-7
Formula:	C ₂₀ H ₁₇ ClN ₂ O ₄
Molecular Weight:	384.82
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	Rebamipideimpurity 4 is an impurity associated with Rebamipide.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5986 mL	12.9931 mL	25.9862 mL
5 mM	0.5197 mL	2.5986 mL	5.1972 mL
10 mM	0.2599 mL	1.2993 mL	2.5986 mL
50 mM	0.052 mL	0.2599 mL	0.5197 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481