

Chick CB11 of type II collagen

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Chick CB11 of type II collagen (Chick cyanogen bromide digested peptide 11) is a cyanogen bromide-digested peptide fragment of type II collagen (CB-peptide) that is useful for studying B cell and T cell epitopes of type II collagen. This product is free of pepsin, making it ideal for investigating T cell epitope specificity through proliferation assays. Additionally, Chick CB11 of type II collagen can serve as an antigen for analyzing antibody epitope specificity via ELISA.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481