

## Dabigatran impurity 10

## Chemical Properties

CAS No. :	91192-27-7
Formula:	C <sub>9</sub> H <sub>8</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight:	176.18
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Dabigatran impurity 75 is a type of Dabigatran impurity.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.676 mL	28.3801 mL	56.7601 mL
5 mM	1.1352 mL	5.676 mL	11.352 mL
10 mM	0.5676 mL	2.838 mL	5.676 mL
50 mM	0.1135 mL	0.5676 mL	1.1352 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481